


M.D.F. EXPLANATORY ONLY

TEST DESK FIG. 1

R.S.R. FIG. 2

TOLL BOARD EXPLANATORY ONLY



N.Z.P.O.		DRAWN		CHECKED		WIRING		CIRCUIT	
AMENDMENT PARTICULARS		DATE		APP.		DATE		APP.	
CANCELED (2222)		7-11-72		B.S.T.		7-11-72		B.S.T.	
NOTE 3 REFERRED TO K CONTACT INSULATION (659)		14-1-54		L.M.K.		11-7-53		D.W.R.	
RECTIFIERS AND E-F RESISTORS HAVE BEEN NOT SHOWN (657)		11-7-53		D.W.R.		11-7-53		D.W.R.	
HOSPITAL CIRCUIT FOR RURAL AUTOMATIC LINES HAMILTON EXCHANGE 2000 TYPE		REF		SIZE		SOV		S	
GBW 10550		REF		SIZE		SOV		S	

- NOTES**
1. FUSE 1-5A PER CIRCUIT
 2. U24 IS STRAPPED TO U11
 3. Ø FOR WIRING SEE GBW 10640 OR EQUIVALENT, FOR R.S.R. & G.B.W. 10870 OR EQUIV. FOR M.D.F.
 4. ALL WIRING NOT OTHERWISE SPECIFIED TO BE 23 SWG. RELAY SET WIRING TO BE 25 S.W.G.
 5. U14 AND U16 TO MAKE CONTACT WHEN THE RELAY SET IS REMOVED.

ETH. BAR 20 S.W.G. R (A) U11
 BATT. FUSE 20 S.W.G. BK (B) U12
 ETH. R U11 TRA9-CO21-COα-CO2-D21-RI-LR2-1α-CR2-T1
 BATT. W U12 RRE-CRe-Ce-YA2-YB2-YC2
 YE2

